IFU


(Instruction for use)

Chave de Mão / Hand Driver / Llave de Mano

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The Medens hand-driver (CMKE-22106) has cylindrical shape and removable anatomical cable, which favors the user's grip. At its active part it has a quick coupling for the fitting of the bone compactor and bone expander tips. The raw material used is stainless steel surgical series 300, according to ASTM F899 303 and aluminum.

APPLICATION

The hand-driver is a surgical instrument indicated to help the user in dental procedures. This product is used manually in conjunction with the bone compactor and bone expander tips, in the preparation of the alveolus for implant installation.

CONTRAINDICATIONS

The hand-driver do not present contraindications, provided they are used correctly for the purposes indicated.

HANDLING

Once sterilized, the hand-driver should be used in compliance with all biosafety and surgical asepsis requirements.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the Medens indications for each instrument.

To attach the tip instrument, hold the hand-driver in the palm of your hand and pull the movable socket sleeve from the end of the quick coupling toward the handle with your fingertips. With the other hand, attach the compactor or expander tip. Release the movable socket sleeve and make sure that the tip is attached to the hand-driver. Take the assembly to the procedure execution point and rotate the cable in the desired direction.

SANITATION AND STERILIZATION

The hand-driver is reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

[Surgical Kit Cleaning Manual](#) 

PRESENTATION

The hand-driver is individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process described above must be followed as soon as it's first used.

PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Before each procedure, make sure that the parts fit together correctly.
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Make sure you have all the necessary instruments for the surgery according to surgical
- planning;
- The use of excessive force can compromise the surgical act;
- Before each procedure, check the conditions of the Medens instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear;
- Always use the sequence of Medens products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Medens from offering any product warranty.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which

can cause damage to the surgical sequence of the implant to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

The Medens hand-driver should be disposed when it loses its functionality.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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