# IFU


(Instruction for use)

## Guia de Brocas / Drill Guide / Guia de Fresas

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological beings and therefore, to obtain results, it must have been improved to the original standards according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

**DESCRIPTION**

The Guide Drill Guides have narrow, regular sizes. Narrows are in lilac and regular in yellow. The internal diameter serves to guide and guide the Guide drills during the osteotomy in the surgical act. The outer diameter serves to be attached to the washer of the surgical guide directing the drill guide according to surgical planning.

The Guide Guide Drills feature a stainless steel surgical cable and guide bush of anodized titanium alloy drills, conforming to ASTM F136. On the rod, there is an indicative recording of the drill bit to be used for the preparation sequence Guide, as shown in the table below:

Product	Drill	Drill Reference	Size	Color
GG-2020	Spear Ø2,0	BLG-2021	Narrow	Lilac
	Spear Ø2,0	BLG-2021	Regular	Yellow
GEG-2832	Ø2,8	BCG-2528 or BLG-2528	Narrow	Lilacs
	Ø3,2	BCG-2832 or BLG-2832	Narrow	Lilac
GRG-2832	Ø2,8	BCG-2528 or BLG-2528	Regular	Yellow
	Ø3,2	BCG-2832 or BLG-2832	Regular	Yellow
GRG-3538	Ø3,5	BCG-3235 or BLG-3235	Regular	Yellow
	Ø3,8	BCG-3538 or BLG-3538	Regular	Yellow

**APPLICATION**

The Guide surgical technique is indicated for surgical procedures in maxilla or mandible bones with a guided instrument sequence, enabling flapless surgery to implant implants.

The Guide Drill Guide is a surgical instrument that serves as a positioning device for Guide drills in regular or narrow washers in the surgical guide prototyped during bone drilling. Compatible with Medens Guide Guided Surgery Technique

**CONTRAINDICATIONS**

These products do not present contraindications provided they are used correctly for the indicated purposes.

**MAGNETIC RESONANCE IMAGE (RM) - INFORMATION**

These products are made of metallic material that can be affected by the RM energy.

For more information, see "Safety Information | Magnetic Resonance Imaging (MRI) "at <https://www.em fi ls.com/ifu>

**HANDLING**

Once sterilized, the Drill Guide Guide should be used to meet all biosafety and surgical asepsis requirements.

All procedures must be carried out by qualified professionals and the handling of the product is inherent to the training of the professional who will be responsible for using the instruments in the best possible way, respecting the indications of the Emphasis for each instrument.

After fixing the Surgical Guide, select the Medens Guide Drill Guide corresponding to the Drill Guide that will be used and according to the washer model. Each Model Drill Guide has a specific Guide Drill Guide,


making sure to use the corresponding assembly correctly. Snap the end of Guide Guide Drill into the washer of the Surgical Guide, holding firmly by the handle. Fit the Guide Drill into the contra-angle and set up the surgical motor with drilling speed as indicated for the selected implant. Take the drill bit to the hole in the Drill Guide. Start the engine and perform the drilling with continuous movements of insertion and removal, with abundant irrigation. This irrigation can be either manual or combined with engine irrigation. Drilling should respect bone quantity and quality, the implant selected, and prior virtual planning. The drilling should be according to the length of the implant and the stop contact of the Drill Guide in the Drill Guide Guide.

Important points in bone perforation:

- To minimize surgical trauma, which is critical to the success of osseointegration of the implant, the sequence of drills should be respected as well as irrigation with physiological serum for bone cooling. Refer to the sequence of the drill bits in the updated In fi ls catalog or at [www.em fi ls.com](http://www.em fi ls.com).
- Avoid interrupting motor rotation with the drill bit inside the surgical cavity.

**SANITATION AND STERILIZATION**

Guide Drill Guides are reusable, supplied non-sterile and must be sanitized and sterilized before and after use, following the recommendations in the document available at the link below:

**PRESENTATION**

Guide Drill Guides are individually packaged in tyvek® surgical paper packaging and clear film display. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

**PRECAUTIONS**

- Do not use if the packaging has been tampered with;
- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the drill guides in their original packaging;
- Be sure to use the drill guide compatible with the drill sequence indicated according to the dimensions of the planned implant;
- Due to the reduced opening of the mouth in the posterior region, the use of surgical drills may be difficult;
- Before each procedure check the conditions of the instruments in fi ls, always respecting their useful life. Replace instruments in case of damage, missing markings, deformation and wear;
- Before each procedure, ensure the perfect fit of the drill guide in the washers of the prototype surgical guide;
- The use of excessive force can compromise the surgical act;
- Always use the product sequence In fi ls. The use of instruments and / or components from other manufacturers does not guarantee the perfect function of the Imphal implant system and disclaims any warranty of the product.

**ADVERSE EFFECTS**

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can lead to unsuccessful perforation and even damage to the patient.

**STORAGE CONDITIONS**

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

**DISPOSAL OF MATERIAL**

The operator must always evaluate the condition of the instrumentation after each use and it must be discarded in accordance with current regulations, complying with the environmental requirements and the criterion of potentially dangerous products, and if the instrument is damaged, markings erased, deformed and worn. , taking into account the different levels of contamination, physical risks and selective collection.

**LIFE CYCLE**

Guide Drill Guides should be discarded when they lose their functionality or the user detects excessive wear.

**SHELF LIFE**

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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