## IFU

(Instruction for use)

## O´ring

PORTUGUÊS

ENGLISH

ESPAÑOL

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**PRAZO DE VALIDADE**

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This product has been designed to integrate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components according to Medens instructions and recommendations. All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

**DESCRIPTION**


The Orings Medens Rings are retention rings that are necessarily connected to other devices, being prosthetic or instrumental component. They are manufactured in high quality autoclavable elastomer. They are non-toxic and amorphous, that is, it is innocuous in the oral cavity.

The Oring Medens Rings are presented in the table below:

CODE	APPLICATION	MATERIAL	COLOR	MAXIMUM WORKING TEMPERATURE(°C)
BMOR-6497	✓ MINI FEMALE O´RING MFOR-4140	POLYURETHANE	BLUE	90°C
AOS-2013	✓ INSTALLATION DRIVERS IHC-XXXX / ILCA-XXXX IQC-XXXX / IQLCA-XXXX ✓ EXPANDERS WITH THREAD EXPR-XXXX ✓ TORQUE WRENCH GEAR ECB-0480 ✓ IMPLANT REMOVER SIE-2030	SILICONE	TRANSLUCENT	200°C
AOS-2108	✓ DRILL EXTENDER PB-2905	SILICONE	TRANSLUCENT	200°C

**APPLICATION**

The BMOR-6497 Ring is mounted together with the MINI FEMALE O´RING (MFOR-4140) component, which in turn fits into the ball-type socket of the retaining abutment and has a retention force necessary for one to support and retain the prosthesis.

The BOR-6004 Ring is mounted together with the key installation tool, thread expander, Torquimetro gear and implant lift. In these instruments, the ring BOR-6004 is responsible for the retention of the square grooves 3.9mm and in the case of the gear, the retention of the lock against angle.

The ring BPB-2108 is used only in the drill extension (PB-2905) and is responsible for retaining the lock against angle.

**CONTRAINDICATIONS**

Retention rings can undergo chemical and physical reactions depending on the environment in which they will be exposed. The table below describes the main characteristics of the raw materials used in the rings and the possible reactions generated in certain situations: See table A.

**HANDLING**

Ask the patient about the need to carry out a professional follow-up and follow the guidelines on care, hygiene and prescription of medications. These guidelines are the responsibility of the professional.

BMOR-6497 - With the frequent removal and insertion of the prosthesis by the patient or even by the lack of parallelism between the fittings, the ring (BMOR-6497) can lose its elastic characteristic or even if it ruptures naturally due to material wear. If this ring needs replacement, remove the prosthesis from the patient’s mouth and with the aid of a curette, completely remove the elastic ring from the inside of the mini female if it is still there. Sanitize the prosthesis before installing the new elastic ring. Install the new elastic ring by pressing it and positioning it with a curette until it is properly seated in the inner channel of the miniature female.

BOR-6004 - With frequent removal and insertion into the fittings, the ring (BOR-6004) can lose its elastic characteristic or even break naturally due to material wear. If this ring is to be replaced with a curette, completely remove the elastic ring from the instrument channel if it is still there. Install the new snap ring by pressing and setting it with a curette until it snaps into place.

BPB-2108 - With frequent removal and insertion into the fittings, the ring (BPB-2108) may lose its elastic characteristic or even break naturally due to material wear. If this ring is to be replaced with a curette, completely remove the elastic ring from the instrument channel if it is still there. Install the new snap ring by pressing and setting it with a curette until it snaps into place.

**FORM OF PRESENTATION AND STERILIZATION**

The rings are individually packaged in tyvek® surgical paper packaging and transparent film display. They are supplied clean because they undergo a process of decontamination, which eliminates any type of dirt, however, they need sterilization before the clinical use.

**PRECAUÇÕES**

- Do not use if the packaging has been tampered with;

- This device must be replaced regularly, every 8 months or when you lose your function.
- Ask the patient to attend the office regularly to check the retention of the prosthesis.
- This product should not be sterilized.
- This product is for single use and can not be reused.

Reusing this product may cause:

Adverse biological effects arising from residues of products, micro-organisms and / or substances from previous uses and / or from reprocessing; changes in the physical, mechanical and chemical characteristics, macro and micro-structural originals of the product that may impair its intended functionality.

- The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of the related products.
- Do not use the product if the packaging has been tampered with.
- Before each procedure, make sure that the parts fit together perfectly.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use original Medens products. The use of prosthetic instruments and / or components from other manufacturers does not guarantee the perfect function of the System and disclaims any warranty of the product.

**ADVERSE EFFECTS**

Adverse effects will only occur if the recommendations for use described are not adhered to.

**STORAGE CONDITIONS**

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage can influence the original characteristics of the packaging and the product.

**DISPOSAL OF MATERIALS**

The professional must always evaluate the conditions of the instruments after each use, and if it is detected that the instrument is damaged, markings erased, deformations and wear, they should be discarded according to the current norms, obeying the environmental requirements and the criterion of potentially dangerous products , taking into account the different levels of contamination, physical risks and selective collection.

**SHELF LIFE**

This product has a useful life of 8 months in function.

**EXPIRATION DATE**

Indicated on the product label.

**TABLE A**

GENERAL	RAW MATERIAL	VITON	POLYURETHANE
	DESIGNATION (ASTMD1418)	FPM	AU
	PROPERTIES	FLUOROUS ELASTOMERS	POLYURETHANE
	ASSIGNMENT	RESISTANCE TO HEAT AND SOLVENTS	ABRASION RESISTANCE
ACID RESISTANCE	DILUTED	GREAT	FEW
	FOCUSED	GREAT	FEW
RESISTANCE TO SOLVENTS	ALIPHATIC HYDROCARBONS (GASOLINE / NAFTALINE)	GREAT	GREAT
	AROMATIC HYDROCARBONS (TOLUOL / XILOL)	GOOD	FEW
	OXYGEN (CETONAS)	FEW	FEW
	LACA SOLVENTS	FEW	FEW
RESISTANCE TO OTHERS	ANIMAL OIL	GREAT	GREAT
	VEGETABLE OIL	GREAT	GREAT
	ABSORPTION OF WATER	GREAT	GOOD
	OXIDATION	GREAT	GREAT
	OZONE	GREAT	GREAT
	HEAT AGING	GREAT	REGULAR
	AGING OF SOLAR LIGHT	GOOD	GOOD
	FLAME (FIRE)	GREAT	FEW
	LOW TEMPERATURES	REGULAR	GREAT
	ELECTRIC INSULATION	GOOD	REGULAR
PHYSICAL CHARACTERISTICS	SPECIFIC DENSITY (RAW MATERIAL)	1.85	1.06
	MINIMUM TEMPERATURE ° C	- 70 ° C	-65Å ° C
	HIGH TEMPERATURE ° C	121 ° C	90 ° C

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