## IFU


(Instruction for use)

### Saca Pilar / Abutment Remover / Saca Pilar

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

## DESCRIPTION

The abutment remover (PSPM-1620) has a cylindrical shape where at its upper tip there is a head with an internal hexagonal geometry for connection in the 1.2 mm driver, which will allow the application of torque to remove the component. In its median region, there is a M2 thread that engages in the component to be withdrawn. The abutment works only in components with morse taper fittings. The raw material used in all the drivers is stainless steel, according to ASTM F899 303.

## APPLICATION

The abutment remover is a surgical or laboratory instrumentation, indicated for bi-digital use in the withdrawal of components that connect only to morse taper implants or analogues. Usually the instrumental is necessary due to the natural interlock existing in the conical system of the morse taper that prevents the removal of the component of the implant even without the screw. The abutment remover is not indicated on components with other connection types other than cone taper. It can be used both in the mouth and in the laboratory for any type of procedure where it is necessary to remove the morse taper component after the application of force to the long axis in the installation.

## CONTRAINDICATIONS

The abutment remover does not present contraindications, provided it's used correctly for the purposes indicated.

## HANDLING

If used in the mouth, the abutment remover should be sterilized and handled only in a sterile environment, by properly trained professionals and in suitable suits at the time of the procedure for withdrawal of the dental morse taper component type. All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the Emfils indications for each instrument.

Collect the abutment using the 1.2 mm hexagonal driver length that best suited to the job to be performed. At the time of use, be sure to align the abutment remover along the long axis of the morse taper component. Bi-digitally, thread the abutment remover in the component clockwise until it becomes detached from the morse taper connection, this may be noticeable because the user notices a small resistance in the threading, indicating that the abutment remover pushed the component out of the connection. By withdrawing the abutment remover, the component is threaded close to it. Loosen the component by unscrewing it counterclockwise.

## SANITATION

The circular scalpels are reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

[Novo Colosso's Surgical Kit Cleaning Manual](#)

## PRESENTATION

The abutment remover is individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as it's first used. Need to be sterilized before using.

## STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

## PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- The use of the abutment remover is indicated only for morse taper components;
- Be aware in cases of patients presenting signs of allergy or hypersensitivity to the chemical components of the instrument material: stainless steel;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, thread loosening or fracture of the prosthetic screw;
- Make sure you have all the necessary instruments to perform the surgery according to surgical planning;
- Before each procedure, check the conditions of

the Emfils instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear;

- Before each procedure, ensure that the abutment remover is in driver correct position, otherwise deformation, scaling or even breaking of the driver can occur in the application of the torque;
- The use of excessive force can compromise the surgical act, damage or lock the driver or component. The abutment remover loose the component very easily, so it should be used only bi-digital force;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty.

## ADVERSE EFFECTS

Adverse effects will only occur if the described use recommendations are not respected.

## STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

## DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

## LIFE CYCLE

The abutment remover should be disposed when it loses its functionality.

## SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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