IFU

(Instruction for use)

Brocas Guide / Guide Drills / Fresas Guide

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological beings and therefore, to obtain results, it must have been improved to the original standards according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.


DESCRIPTION

The Guide surgical drill bits that make up the Medens Guide surgical kit form a set of 10 drill bits to drill the bone tissue for implant installation and / or surgical guide fixation pin. They are available in 3 models: Ø1,3mm pin, drill bit and helical drills (short and long).

The drill bit of the fixing pin has an active end with a diameter of 1.3mm and at the other end has a counter-angle fitting.

The drill launches at its active end, has spear-shaped tip with high cutting power. At the other end the boom drill has a counter-angle fitting.

The helical drills are staggered and have at their active end a cylindrical profile with cutting geometry compatible with the macrostructure of the implants. At the tip of the active end we have a smaller cylindrical geometry, with the same diameter of the previous drill, which serves as a guide for the posterior drilling. At the other end the helical drills also have the counter-angle fitting. The laser markings take into account the actual implant length, with an additional 11 mm length to compensate for the soft tissue thickness and height of the washer and the Drill Guide. See image below:

| Nominal mark height (mm) | Actual mark height (mm) |
|--------------------------|-------------------------|
| 7 | 18 |
| 8,5 | 19,5 |
| 10 | 21 |
| 11,5 | 22,5 |
| 13 | 24 |
| 14,5 | 25,5 |
| 16 | 27 |

Guide surgical bits receive the application of a thin biocompatible carbon coating that imparts dark coloration to the cutting surface. This coating has as main purposes:

- Increased resistance to oxidation;
- Reduced friction between drill and bone;

- Increased wear resistance;
- Reduction of the heat generated in the bone during the osteotomy;
- Optimization of the visualization of height markings;

The Guide surgical drills that make up the Medens surgical kit are presented in the table below:

| Sales Code | Model | Drilling Diameter (mm) |
|------------|---------------------------------|------------------------|
| BFG-1335 | DRILLING PIN DRIVE Ø1,3MM GUIDE | 1,3 |
| BLG-2021 | BROCA LAUNCHES Ø2,0MM GUIDE | 2,0 |
| BCG-2528 | SHORT DRIVE Ø2,5 / 2,8MM GUIDE | 2,8 |
| BCG-2832 | SHORT DRIVE Ø2,8 / 3,2MM GUIDE | 3,2 |
| BCG-3235 | SHORT DRIVE Ø3,2 / 3,5MM GUIDE | 3,5 |
| BCG-3538 | SHORT DRIVE Ø3,5 / 3,8MM GUIDE | 3,8 |
| BLG-2528 | LONG DRILL Ø2,5 / 2,8MM GUIDE | 2,8 |
| BLG-2832 | LONG DRILL Ø2,8 / 3,2MM GUIDE | 3,2 |
| BLG-3235 | LONG BRUSH Ø3,2 / 3,5MM GUIDE | 3,5 |
| BLG-3538 | LONG DRILL Ø3,5 / 3,8MM GUIDE | 3,8 |

The raw material used for the Guide surgical drills is tempered stainless steel surgical steel, according to ASTM F899.

APPLICATION

The Medens Guide surgical technique is indicated for surgical procedures in maxilla or mandible bones with a guided instrument sequence, enabling flapless surgery to implant implants in edentulous areas.

The guide drill 1.3 is an instrument used for bone perforation allowing the insertion of the pin of fixation in the Prototype Guide. The Fixing Pin serves as a fixation bracket for the Surgical Guide during the implant installation procedure.

The Guide spear drill is an instrumental used for the demarcation and the rupture of the cortical bone.

Guide drills are used for bone drilling in the preparation of the artificial surgical alveolus for implant installation. At the end it will result in a staggered drilling that will have a surgical store suitable for the macro geometry of the implants.

CONTRAINDICATIONS

Guide surgical drills do not have contraindications, provided they are used correctly for the indicated purposes.

MAGNETIC RESONANCE IMAGE (RM) - INFORMATION

These products are made of metallic material that can be affected by the RM energy.

For more information, see "Safety Information | Magnetic Resonance Imaging (MRI)" at <https://www.Medens.com/ifu>

HANDLING

Once sterilized, the Guide drill should be used to meet all biosafety and surgical asepsis requirements.

All procedures must be carried out by qualified professionals and the handling of the product is inherent to the training of the professional who will be responsible for using the instruments in the best possible way, respecting the indications of the Emphasis for each instrument.

For drilling with Broca Guide, after fixing the Surgical Guide in the oral cavity, select the Medens Guide Drill Guide corresponding to the Drill Guide that will be used and according to the washer model. Each Broca Guide has a specific Medens Guide Drill Guide. Snap the end of the Drill Guide into the washer of the Surgical Guide, holding the cable firmly. Fit the Guide Drill into the contra-angle and set up the surgical motor with drilling speed as indicated for the selected implant. Take the drill bit to the hole in the Drill Guide. Start the engine and drill with continuous insertion and removal movements with abundant irrigation. This irrigation can be either manual or combined with engine irrigation. Drilling should respect bone quantity and quality, the implant selected, and prior virtual planning. The drilling should be according to the length of the implant and the laser markings of the drill bit. The laser markings on the drill bit represent the drill depth for the length of the implant.

Except drill 1.3, other drills should be used in conjunction with drill guide and surgical guide. The fixing pin drill will be used directly on the prototyped surgical guide or the washer, ensuring greater accuracy.

Important points in bone perforation:

To minimize surgical trauma, which is critical to the success of osseointegration of the implant, the sequence of drills should be respected as well as irrigation with physiological serum for bone cooling. Refer to the sequence of the drill bits in the updated In fi ls catalog or at www.em fi ls.com.

Avoid interrupting motor rotation with the drill bit inside the surgical cavity.

SANITATION AND STERILIZATION

Guide surgical drills are reusable, supplied non-sterile and must be sanitized and sterilized before and after use, following the recommendations in the document available at the link below:

[Surgical Kit Cleaning Manual](#) 

PRESENTATION

Guide surgical drills are individually packaged in tyvek® surgical paper packaging and clear film display. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

PRECAUTIONS

- Do not use if the packaging has been tampered with;
- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the drills in their original packaging;
- Care must be taken in the hygiene of the drills because it is a sharp metal instrument, so as not to suffer any type of physical trauma;
- Be sure to use the drill bit compatible with the specified bit sequence according to the dimensions of the planned implant;
- The drills should not be re-drilled;
- Failure to replace the drill, as recommended by Emphil, may lead to undue bone heating, compromising the success of the procedure. The drills should not be compromised. Refer to the "Lifetime" item at the end of this instruction manual;
- Due to the reduced opening of the mouth in the posterior region, the use of surgical drills may be difficult;
- Before each procedure check the conditions of the instruments in fi ls, always respecting their useful life. Replace instruments in case of

damage, missing markings, deformation and wear;

- Before each procedure, make sure that the drill is perfectly in place against the angle, otherwise deformation or even breaking of the drill bit may occur;
- The use of excessive force can compromise the surgical act;
- Always use the product sequence In fi ls. The use of instruments and / or components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any product warranty

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can lead to unsuccessful perforation and even damage to the patient.

STORAGE CONDITIONS

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use, and if the instrument is detected to be damaged, markings erased, deformations and wear, should be discarded.

LIFE CYCLE

This product is recommended for up to 30 uses, subject to the conditions of use recommended by Emfs. Regardless of the number of times the instrument was used, the operator should always evaluate the instrument's conditions after each use.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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