IFU


(Instruction for use)

Montador / Assembler / Ensamblador

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

DESCRIPTION

The EMFILS assembler (MMPA-4800) facilitates the handling and installation of the EMFILS mini angled abutment by transporting the package to the installed implant or laboratory analogue. When used during the surgical procedure, it minimizes the risk of contamination as there is less manual contact with the product and, moreover, assists in the adaptation, alignment and initial retention of the angled mini pillar. The raw material used in the EMFILS assembler is stainless steel surgical steel according to ASTM F899.

APPLICATION

The EMFILS assembler (MMPA-4800) is indicated for transport and installation of the mini-column abutment in the installed implant (surgical procedure) or laboratory analogue (laboratory procedure). It can also be used as an alignment reference for the insertion angle. On one end it has a 1.4 thread to connect with the thread of the standard mini abutment.

It should not be used on straight mini pillars, only on angled ones.

CONTRAINDICATIONS

The EMFILS assembler (MPA-4800) does not present any contraindications, provided it is used correctly for the indicated purposes.

HANDLING

If used in the surgical procedure, the EMFILS assembler (MMPA-4800) should be sterilized before use. Once sterilized, it should be handled only in a sterile environment, by properly trained professionals and in appropriate attire at the time of surgery for the installation of dental implants.

All procedures must be carried out by qualified professionals and the handling of the product is inherent in the training of the professional, who will be responsible for using the instruments in the best possible way, respecting the indications of the Emphasis for each instrument.

Important points in use:

- Avoid side loading.
- The product has no torque control.
- Excessive force can damage the product.

SANITATION AND STERILIZATION

The EMFILS assembler (MMPA-4800) is reusable, supplied non-sterile and should be sanitized and sterilized before and after use, when used in a surgical procedure, following the recommendations of the document available at the link below:

[Novo Colosso's Surgical Kit Cleaning Manual](#)

PRESENTATION

The product is individually packaged in tyvek® surgical paper packaging and clear film display. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

PRECAUTIONS

- Do not use if the packaging has been tampered with;
- When used in a surgical procedure, this product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Avoid side load on the product;
- The product has no torque control.
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Before each procedure check the conditions of the instruments In fi ls, always respecting their useful life. Replace instruments in case of damage, missing markings, deformation and wear;
- The use of excessive force can compromise the surgical act and damage the product or component connected to it;
- Always use the product sequence In fi ls. The use of instrumentation and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any warranty of the product.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instruments is inadequate, which can cause damage to the component to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

DISPOSAL OF MATERIAL

The operator must always evaluate the condition of the instrumentation after each use and it must be discarded in accordance with current regulations, complying with the environmental requirements and the criterion of potentially dangerous products, and if the instrument is damaged, markings erased, deformed and worn, taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

The product should be discarded when it loses its functionality or the wearer detects excessive wear.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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