IFU


(Instruction for use)

Sistema CADCAM / CADCAM System / Sistema CADCAM

PORTUGUÊS

ENGLISH

ESPAÑOL

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This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The CAD-CAM system for digital dentistry allows customization of the prostheses in drawing through a specific software and the subsequent machining of these prosthesis in a milling machine connected in a computer. To use this system, the user must request the Emfils library corresponding to the software platform that he has. This library contains all types of interfaces in the Emfils CAD-CAM system and meets to any kind of work that the user prefers to perform. Download your corresponding library on Emfils' own website or request it from your seller. The Emfils

CAD-CAM system has the following physical components:

Scan Body or Jig:

Of transitory use, the scan body or jig has different interfaces depending on the type of implant or base abutment chosen in the treatment. It works as a transfer to the virtual environment that identifies the exact three-dimensional position of the implant installed in the mouth, through the trapezoidal geometry in its superior part, which creates reference faces that will guide you after the scan. It can be used both in the plaster model (flatbed scanner) on the laboratory, and in the mouth (intraoral scanner) on the dental clinic. The fixing of the scan body in the component must be by a specific clinical screw.

The raw material of the scan body is the titanium alloy (grade 5), according to ASTM F136. This titanium alloy is blasted that makes the surface less shiny, which favors the precision of the scanning.

Ti-Base or titanium base:

Of definitive use, the ti-base or titanium base has different interfaces depending on the type of implant or base abutment chosen in the treatment. It works as an intermediate component between the implant or base abutment and the prosthesis. The ti-base or titanium base fixes the prosthesis through dental resin cement. It has a cylindrical geometry with variations in diameter, height (gingival heights) and can have a non-rotational interface for single prosthesis or rotational interface for multiple prosthesis, depending on the diameter and model of the installed component. The height of the cementation area is constant for all ti-base interfaces: 4 mm and also have external retention (groove and milling), which prevents both loosening and turning of the prosthesis.

The raw material of the ti-base is the titanium alloy (grade 5), according to ASTM F136.

The scan bodies of the CAD-CAM system are presented in the table below:

Sales Code	Interface Type								Implant or Base Abutment Diameter	Interface		Screw
	MT	ZL	HE	IH	Base Abutment					Rotat.	Non Rotat	
					Multi	Mini 4,0-5,0	Mini 4,8	Angled				
JCM-3512	X								3,5-4,0-5,0		X	M1.6
JHE-3512			X						3,5-4,0-5,0		X	M2
JHI-3512				X					3,3-3,5-4,0-5,0		X	M2
JLZ-3512		X							3,5-4,0-5,0		X	M1.6
JMP-4012						X			4,0-5,0	X		M1.6
JMPP-4812							X	X	4,8	X		M1.4
JNRMU-4012					X				4,0-5,0		X	M2

The ti-bases of the CAD-CAM system are presented in the table below:

Sales Code	Interface Type								Implant or Base Abutment Diameter	Height	Interface		Screw
	MT	ZL	EH	IH	Base Abutment						Rotat.	Non Rotat	
					Multi	Mini 4,0-5,0	Mini 4,8	Angled					
LKLMPP-4805							X	X	4,8	---	X		M1.4
LKLNNR-4005					X				4,0-5,0	---		X	M2
LKLRMP-4005						X			4,0-5,0	---	X		M1.6
LKLRNR-4005					X				4,0-5,0	---	X		M2
LKNLE-3504			X						3,5	---		X	M2
LKNLE-4004			X						4,0-5,0	---		X	M2
LKNLI-3504		X							3,5	---		X	M1.6
LKNLI-4004		X							4,0-5,0	---		X	M1.6
LKNLSI-3504				X					3,3-3,5	---		X	M2
LKNLSI-4004				X					4,0-5,0	---		X	M2
LKRLE-3504			X						3,5	---	X		M2
LKRLE-4004			X						4,0-5,0	---	X		M2
LKRLI-3504		X							3,5	---	X		M1.6
LKRLI-4004		X							4,0-5,0	---	X		M1.6
LKRLSI-3504				X					3,3-3,5	---	X		M2
LKRLSI-4004				X					4,0-5,0	---	X		M2
LKRNL-4015	X								3,5-4,0-5,0	1,5		X	M1.6
LKRNL-4030	X								3,5-4,0-5,0	3,0		X	M1.6
LKRNL-4045	X								3,5-4,0-5,0	4,5		X	M1.6

The Emfils may not have developed the library for some specific software platforms on the market, so before installing or requesting, check which platform your software works on.

APPLICATION

Scan Body or Jig:

This prosthetic component is used to transfer the three-dimensional position of the implant or base abutment installed in the mouth. This transfer can be from the analogue, during the scanning of a plaster model through a flatbed scanner, or the scanning of the installed component itself in the mouth through an intraoral scanner in the clinic. The user must always plan each case to use the CAD-CAM technique, respecting the tolerances of the components, in view of the aesthetics and functionality. Before finalizing the work, make sure that the screw has been tightened to the required torque.

Ti-Base or titanium base:

This prosthetic component is indicated according to the implant or base abutment model. On it, can be cemented a customized abutment for a cemented crown, or can be cemented a finished crown for attachment through the screw. The user must always plan each case to use the CAD-CAM technique, respecting the tolerances of the components, in view of the aesthetics and functionality. Before finalizing the work, make sure that the screw has been tightened to the required torque.

CONTRAINDICATIONS

In the case of ti-bases, reuse is not recommended and may directly bad influences to the quality and success of the procedure.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.emfils.com.br/ifu>

HANDLING

Scan Body or Jig:

In case of the plaster model scanning, install it on the analogue of the same connection with the corresponding clinical screw, using a 1.2 mm hexagonal driver and applying a maximum bi-digital torque of 10 N.cm. Check if the component is properly seated in the analogue because there must be no rotational or vertical clearance. Clearances indicate that the scanning body is not properly seated and should be repositioned. Take the set to the flatbed scanner used for the procedure

following the instructions of the equipment supplier. If you use a developer spray for a better scan preview, avoid excessive spraying on the scan body because it has a matte surface prepared for the process.

In the case of direct intraoral scanning, install the scan body in the mouth on the implant or base abutment of same connection with the corresponding clinical screw, using a 1.2 mm hexagonal driver and applying a maximum bi-digital torque of 10 N.cm. Check if the component is properly seated in the implant or base abutment because there must be no rotational or vertical clearance. Clearances indicate that the scanning body is not properly seated and should be repositioned. Start scanning with the intraoral scanner, following the instructions of the equipment supplier.

Ti-Base or titanium base:

Any ti-base of the line does not allow customization of the cementing area because they already have a maximum height of 4 mm and must remain intact. Make the crown or structure using your CAD software with the compatible platform to work together with the Emfils digital library that consist of 3D set files containing the milling matrix required for the interface geometry used in the job.

Position the ti-base over the plaster model and tighten the set screw over it with a slight tightening using the 1.2 mm hex driver. Protect the access of the screw during the process of cementation of the crown or structure. The handling of the cement should follow the manufacturer’s cement instructions. Apply the cement on the ti-base outer portion and

press the piece into place by fitting it perfectly on the ti-base. The excess cement that is spilled out of the hole should be removed immediately. After the cement is secured, the infrastructure is unscrewed from the analog and the excess cement remaining at the edge of ti-base must be removed. It is recommended to use chemically activated resin cements for metal bonding. The ti-base cementation on the work is the conclusive step of the process.

PRESENTATION

All parts that make up the CAD-CAM system are individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Both the scanning body and the ti-base are supplied non-sterile. All ti-bases models are indicated for single use only. The scan bodies are for single use only if scanned in the mouth (intraoral scanner). If they are scanned on the plaster model (flatbed scanner), can be used more than once, as long as the connection isn’t with excessive wear or damaged, so that to compromise their positioning when scanning.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Before installing the products, make sure that they have the appropriate form of use for the selected prosthetic system;
- Special care should be taken during prosthesis design for what the interocclusal space is not insufficient or excessive.
- Ensure that the parts are not swallowed or aspirated by the person who is receiving the product (patient);
- The scan bodies are for single use only if scanned in the mouth (intraoral scanner). If they are scanned on the plaster model (flatbed scanner), they can be used more than once, as long as the connection isn’t with excessive wear or damaged, so that to compromise their positioning when scanning.
- Using the scanner body with the damaged connection, platform, or geometry may lead to unsatisfactory scanning results;
- All ti-bases models are indicated for single use only. Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can

cause damage to the material and make the system unusable;

- Carefully verify that the scan body is well-fitting to the model, implant or base abutment analog, without rotational or vertical clearances;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- When installing components on the implant, base abutment or analogue, be sure to align the connection and/or the screw, thus avoiding damage inside the piece or locking the system;
- The ti-base is not indicated for mouth cementation;
- The ti-base for CAD-CAM system must be perfectly seated to the implant or base abutment after the necessary torque has been applied to the screw. This check is recommended using periapical radiographs;

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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