## IFU


(Instruction for use)

### Pilar Provisório / Temporary Abutment / Pilar Provisional

PORTUGUÊS

ENGLISH

ESPAÑOL

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Código de Venda	Tipo Conexão								Tipo Prótese		Diâmetro Implante ou Transmuc.	Altura Cinta	Altura Cimentação	Parafuso
	CM	LZ	HE	HI	Transmucoso				Unit.	Múlt.				
					Multi	Mini 4,0-5,0	Mini 4,8	Angul. 4,8						
CPNE-3511			X						X		3,5	---	9,0	M2
CPNE-4011			X						X		4,0	---	9,0	M2
CPNE-5011			X						X		5,0	---	9,0	M2
CPNI-3511		X							X		3,5	---	9,0	M1.6
CPNI-4011		X							X		4,0	---	9,0	M1.6
CPNI-5011		X							X		5,0	---	9,0	M1.6
CPNNR-4011					X				X		4,0	---	10,0	M2
CPNNR-5011					X				X		5,0	---	10,0	M2
CPPNR-3314				X					X		3,3-3,5	---	15,0	M2
CPPNR-4014				X					X		4,0	---	15,0	M2
CPPNR-5014				X					X		5,0	---	15,0	M2
CPPR-3314				X						X	3,3-3,5	---	15,0	M2
CPPR-4014				X						X	4,0	---	15,0	M2
CPPR-5014				X						X	5,0	---	15,0	M2
CPRE-3511			X							X	3,5	---	9,0	M2
CPRE-4011			X							X	4,0	---	9,0	M2
CPRE-5011			X							X	5,0	---	9,0	M2
CPRI-3511		X								X	3,5	---	9,0	M1.6
CPRI-4011		X								X	4,0	---	9,0	M1.6
CPRI-5011		X								X	5,0	---	9,0	M1.6
CPRMP-4011						X				X	4,0	---	10,0	M1.6
CPRMP-5011						X				X	5,0	---	10,0	M1.6
CPRNM-3507	X								X		3,5-4,0-5,0	1,0	9,0	M1.6
CPRNM-3515	X								X		3,5-4,0-5,0	1,5	9,0	M1.6
CPRNM-3530	X								X		3,5-4,0-5,0	3,0	9,0	M1.6
CPRNM-3545	X								X		3,5-4,0-5,0	4,5	9,0	M1.6
CPRNM-3560	X								X		3,5-4,0-5,0	6,0	9,0	M1.6
CPRNR-4011					X					X	4,0	---	10,0	M2
CPRNR-5011					X					X	5,0	---	10,0	M2
PMPP-4810							X	X		X	4,8	---	10,0	M1.4

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

## DESCRIPTION

The temporary abutments are part of the prosthetic components family and are used as an intermediate between the implant or base abutment and the prosthesis (crown) that must be cemented. The product is available with different fittings, depending on the type of implant or base abutment used and is intended for the manufacture of temporary screwed prostheses.

The abutments have a cylindrical surface with retentive grooves and flat surfaces for accomplishment of the waxing or direct cementation with commercial prosthesis, without the concern of work displacement.

The raw material used is the titanium alloy (grade 5), according to ASTM F136.

The temporary abutments are presented in the table below:

Sales Code	Connection Type								Prosthesis		Implant Diameter or Base Abut	Collar Height	Cimentation Height	Screw
	MT	ZL	EH	IH	Base Abutments				Unit.	Múlt.				
					Multi	Mini 4,0-5,0	Mini 4,8	Angled 4,8						
CPNE-3511			X						X		3,5	---	9,0	M2
CPNE-4011			X						X		4,0	---	9,0	M2
CPNE-5011			X						X		5,0	---	9,0	M2
CPNI-3511		X							X		3,5	---	9,0	M1.6
CPNI-4011		X							X		4,0	---	9,0	M1.6
CPNI-5011		X							X		5,0	---	9,0	M1.6
CPNNR-4011					X				X		4,0	---	10,0	M2
CPNNR-5011					X				X		5,0	---	10,0	M2
CPPNR-3314				X					X		3,3-3,5	---	15,0	M2
CPPNR-4014				X					X		4,0	---	15,0	M2
CPPNR-5014				X					X		5,0	---	15,0	M2
CPPR-3314				X						X	3,3-3,5	---	15,0	M2
CPPR-4014				X						X	4,0	---	15,0	M2
CPPR-5014				X						X	5,0	---	15,0	M2
CPRE-3511			X						X		3,5	---	9,0	M2
CPRE-4011			X						X		4,0	---	9,0	M2
CPRE-5011			X						X		5,0	---	9,0	M2
CPRI-3511		X								X	3,5	---	9,0	M1.6
CPRI-4011		X								X	4,0	---	9,0	M1.6
CPRI-5011		X								X	5,0	---	9,0	M1.6
CPRMP-4011						X				X	4,0	---	10,0	M1.6
CPRMP-5011						X				X	5,0	---	10,0	M1.6
CPRNM-3507	X								X		3,5-4,0-5,0	1,0	9,0	M1.6
CPRNM-3515	X								X		3,5-4,0-5,0	1,5	9,0	M1.6
CPRNM-3530	X								X		3,5-4,0-5,0	3,0	9,0	M1.6
CPRNM-3545	X								X		3,5-4,0-5,0	4,5	9,0	M1.6
CPRNM-3560	X								X		3,5-4,0-5,0	6,0	9,0	M1.6
CPRNR-4011					X					X	4,0	---	10,0	M2
CPRNR-5011					X					X	5,0	---	10,0	M2
PMPP-4810							X	X		X	4,8	---	10,0	M1.4

## APPLICATION

The temporary abutments are indicated for temporary jobs in unitary prostheses (non-rotational fitting) or multiple (rotational fitting) screwed onto the implant or base abutment, according to available interocclusal space, existing gingival height and three-dimensional position of the implant. Indicated to be used in immediate rehabilitation procedures, in the maxilla or mandible. While choosing the temporary abutment, make sure how to use the selected component. Before finishing the clinical work, make sure that the screw has

been tightened with the necessary torque.

## CONTRAINDICATIONS

The temporary restoration made on the temporary abutment should not be placed in occlusion.

Contraindicated in cases where the interocclusal space is insufficient and / or the three-dimensional position of the implant is unsatisfactory.

## MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.emfils.com.br/ifu>

## HANDLING

It can be used in two ways:

- a) with the appropriate connection, fit the screw inside the abutment, install it over the implant and prepare it after application of the recommended torque on the screw;
- b) transfer the implant installed to the working model (plaster). With the appropriate connection, fit the screw inside the selected abutment, position it on the corresponding analogue in the model and prepare it.

The temporary abutments should be used in parallel prostheses without angulation. The cementation of the prosthesis is performed in the mouth following the instructions of the cement manufacturer. In this case, it is worth remembering that it is necessary to know the technique to avoid excess cement and chemical burns in the tissues adjacent to the implant or base abutment. Protect screw access during the cementation process, if applicable.

It is recommended to use stainless steel screw for lab work and the titanium alloy screw for fixing in the mouth, always taking into account the recommended torque.

### PRESENTATION

The abutments are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

### STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

### PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Before installing the product, make sure that it has the appropriate the way in which of use for the selected prosthetic component.
- Ensure that the pieces are not swallowed

or aspirated by the person receiving the product (patient);

- Any type of abutment is for single use only and should not be re-sterilized.
- Reusing this product may cause:
  - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
  - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
  - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Special care should be taken during prosthesis design for what the interocclusal space is not insufficient or excessive;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;
- The component must be perfectly seated to the implant, base abutment or analog after the necessary torque has been applied to the screw. This check is recommended using periapical radiographs with the parallelism technique, when applicable.
- When using a unitary restoration, it is medically advisable to maintain a minimum foundation height of 4 mm above the gum;
- It is recommended that the temporary prosthesis be hygienised before installation in the mouth;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

### ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

### STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

### DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

### DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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Código de Venda	Tipo Conexão								Tipo Prótese		Diâmetro Implante ou Transmuc.	Altura Cinta	Altura Cimentação	Parafuso
	CM	LZ	HE	HI	Transmucoso				Unit.	Múlt.				
					Multi	Mini 4,0-5,0	Mini 4,8	Angul. 4,8						
CPNE-3511			X						X		3,5	---	9,0	M2
CPNE-4011			X						X		4,0	---	9,0	M2
CPNE-5011			X						X		5,0	---	9,0	M2
CPNI-3511		X							X		3,5	---	9,0	M1.6
CPNI-4011		X							X		4,0	---	9,0	M1.6
CPNI-5011		X							X		5,0	---	9,0	M1.6
CPNNR-4011					X				X		4,0	---	10,0	M2
CPNNR-5011					X				X		5,0	---	10,0	M2
CPPNR-3314				X					X		3,3-3,5	---	15,0	M2
CPPNR-4014				X					X		4,0	---	15,0	M2
CPPNR-5014				X					X		5,0	---	15,0	M2
CPPR-3314				X						X	3,3-3,5	---	15,0	M2
CPPR-4014				X						X	4,0	---	15,0	M2
CPPR-5014				X						X	5,0	---	15,0	M2
CPRE-3511			X						X		3,5	---	9,0	M2
CPRE-4011			X						X		4,0	---	9,0	M2
CPRE-5011			X						X		5,0	---	9,0	M2
CPRI-3511		X								X	3,5	---	9,0	M1.6
CPRI-4011		X								X	4,0	---	9,0	M1.6
CPRI-5011		X								X	5,0	---	9,0	M1.6
CPRMP-4011						X				X	4,0	---	10,0	M1.6
CPRMP-5011						X				X	5,0	---	10,0	M1.6
CPRNM-3507	X								X		3,5-4,0-5,0	1,0	9,0	M1.6
CPRNM-3515	X								X		3,5-4,0-5,0	1,5	9,0	M1.6
CPRNM-3530	X								X		3,5-4,0-5,0	3,0	9,0	M1.6
CPRNM-3545	X								X		3,5-4,0-5,0	4,5	9,0	M1.6
CPRNM-3560	X								X		3,5-4,0-5,0	6,0	9,0	M1.6
CPRNR-4011					X					X	4,0	---	10,0	M2
CPRNR-5011					X					X	5,0	---	10,0	M2
PMPP-4810							X	X		X	4,8	---	10,0	M1.4

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