# IFU

(Instruction for use)

## Brocas / Drills / Fresas

PORTUGUÊS

ENGLISH

ESPAÑOL

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#### **PRAZO DE VALIDADE**

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
This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

**DESCRIPTION**

The Diamond surgical drill bits that make up the Medens surgical kit form a set of 8 drills designed to pierce the bone tissue to implant the implant. They are presented in 2 models: The drill spears and the helical drills. The drill launches at its active end, has spear-shaped tip with high cutting power. At the other end the boom drill has a counter-angle fitting. The stepped helical drills have at their active end a cylindrical profile and cutting geometry compatible with the macrostructure of the implants. At the tip of the active end we have a smaller cylindrical geometry, with the same diameter of the previous drill, which serves as a guide for the posterior drilling. At the other end the helical drills also have the counter-angle fitting. All drill bits have laser markings that determine the depth of the drilling according to the lengths of the implants (see image below). In its central part, the drills have a larger diameter ring, which serves as the mechanical stop of the depth limiter.

Diamond surgical drill bits receive the application of a thin biocompatible carbon coating that imparts dark coloration to the cutting surface. This film has as main purposes:

- Increased resistance to oxidation;
- Reduced friction between drill and bone;
- Increased wear resistance;
- Reduction of the heat generated in the bone during the osteotomy;
- Optimization of the visualization of height markings;

The Diamond surgical drills that make up the Medens surgical kit are presented in the table below:

Sales Code	Model	Drill Diameter	Surgical Protocol Implant Ø	
			Bone Type I and II	Bone Type III and IV
BRD-2021	Initial Drill	2,0	2,5	---
BRD-2025	Stepped Twist Drill	2,5	3,0	---
BRD-2528	Stepped Twist Drill	2,8	---	3,3-3,5
BRD-2832	Stepped Twist Drill	3,2	3,3-3,5	---
BRD-3235	Stepped Twist Drill	3,5	---	4,0
BRD-3538	Stepped Twist Drill	3,8	4,0	---
BRD-3845	Stepped Twist Drill	4,5	---	5,0
BRD-4548	Stepped Twist Drill	4,8	5,0	---

The raw material used for Diamond surgical drill bits is tempered stainless steel surgical steel, conforming to ASTM F899.

**APPLICATION**

The Diamond initial drill is an instrument used for the demarcation and disruption of the cortical bone. Diamond stepped twist drills are used for bone drilling in the preparation of the artificial surgical alveolus for implant installation. The twist drills also work as pilot drills because they have a smaller cylindrical geometry at the end of the active tip, with the same diameter as the previous drill, which serves as a guide for later drilling. At the end, we will have a staggered drilling that will have a surgical alveolus suitable for the macro geometry of Medens implants.

**CONTRAINDICATIONS**

Diamond surgical drills do not present contraindications, provided they are used correctly for the purposes indicated.

**MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION**

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.medens.com/ifu>

**HANDLING**

Once sterilized, the Diamond surgical drills should be used in compliance with all biosafety and surgical asepsis requirements.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the Medens indications for each instrument.

Insert the drill into the counter-angle fitting and specific the rotation. Start the engine and drill the bone until the desired depth.

Important points in a bone drilling:

- o minimize surgical trauma, which is critical to the success of osseointegration of the implant, the sequence of drills must be respected, as well as irrigation with saline for bone cooling. See the drill sequencing in the updated Medens catalog or at [www.Medens.com](http://www.Medens.com).
- Avoid interrupting motor rotation with the drill bit inside the surgical cavity.

**SANITATION AND STERILIZATION**

Diamond surgical drills are reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

Surgical Kit Cleaning Manual

**PRESENTATION**

Surgical drills are individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as it’s first used.

**PRECAUTIONS**

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;

- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Should be careful in the drill sanitation process because it is a sharp metal instrument, avoiding any type of physical trauma;
- Be sure to use the drill compatible with the specified drill sequence according to the dimensions of the planned implant;
- Surgical drills should not be re-sharped;
- The non-replacement of a drill, as recommended by Medens, can generate undue bone heating, compromising the success of the procedure. The drills should not be with their sharpening committed. Refer to the "Life Cycle" section at the end of this instruction manual;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Before each procedure, check the conditions of the Medens instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear.
- Before each procedure, make sure that the drill is perfectly coupled into the counter-angle, otherwise deformation or even breaking of the drill may occur;
- The use of excessive force can compromise the surgical act;
- Always use the sequence of Medens products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Medens from offering any product warranty.

#### **ADVERSE EFFECTS**

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can lead to unsuccessful of the drilling and even damage to the patient.

#### **STORAGE CONDITIONS**

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

#### **DISPOSAL OF MATERIAL**

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

#### **LIFE CYCLE**

This product is recommended for up to 30 uses, provided that Medens-recommended conditions of use are followed. Regardless of the number of times the instrument has been used, the professional must always evaluate the instrument's conditions after each use.

#### **SHELF LIFE**

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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