IFU

(Instruction for use)

Chave Hexagonal Bidigital / Bidigital Hexagonal Driver / Llave Hexagonal Bidigital

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The Medens bi-digital hexagonal drivers come in 2 lengths, short and long. At its active tip, it has a hexagonal geometry of 1.2 mm. At the other tip, the drivers have an anatomically designed body for bi-digital use, and has an independently-driven washer and hole for dental floss installation, to increase safety during use. Although the washer has circular movements independent of the body, it is fixed and can't be disassembled from the instrument.

The raw material used is tempered stainless steel, according to ASTM F899.

The Medens bi-digital hexagon drivers are presented in the table below:

Código de Venda	Modelo		Comprim. Útil
	Curto	Longo	
CHCC-1216	X		11,0
CHLC-1226		X	20,0

A distinctive feature found in our 1.2 mm hexagonal installation keys is the Collection System, which enables the efficient fastening of the screw at the tip of the key, preventing screw drop, contamination, or even risks of the screw being aspirated or swallowed by the patient.

APPLICATION

The Medens bi-digital drivers are indicated to help in the installation of all screws and some prosthetic components of the Medens system line, that have a 1.2mm internal hexagon. The bi-digital hexagonal driver can be used as surgical or laboratory instrument.

CONTRAINDICATIONS

The bi-digital drivers do not present contraindications, provided they are used correctly for the purposes indicated.

HANDLING

Once sterilized, the bi-digital hexagonal driver should be used in compliance with all biosafety and surgical asepsis requirements.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the instruments in the best possible way, respecting the Medens indications for each instrument.


Before using the product, it is recommended to tie a dental floss into the hole in the washer, to avoid the risk that the patient can swallow or suck the part.

To capture the screw or prosthetic component,

fitting the active tip of the bi-digital driver into the hexagon of the piece. Press the bi-digital driver, make sure the correct fit, and move it clockwise for insertion, applying a slight squeeze. Make counterclockwise movements for removal.

SANITATION AND STERILIZATION

The bi-digital hexagonal drivers are reusable, provided non-sterile and should be sanitized and sterilized before and after use by following the steps mentioned in the link below:

PRESENTATION

The Medens bi-digital hexagon drivers are individually packaged in tyvek® surgical paper and transparent film display. The instrument undergoes a decontamination process that eliminates any type of dirtiness, but the sanitation process, described above, must be followed as soon as it's first used.

PRECAUTIONS

- Do not use if the product packaging is broken;
- This product should be used sterile. Do not sterilize under dry heat (stove). Always use steam and pressure autoclave;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- Avoid side load on the switch;
- The drivers have no torque control;
- Because of reduced mouth opening in the posterior region, the use of long driver associated with prosthetic components in premolar and molar regions is not recommended;
- Maximum care is recommended when using the instruments so that they are not swallowed or aspirated by the patient. Always use the dental floss technique on the washer hole at the top of the scalpel, giving greater safety to the procedure;
- Before each procedure, check the conditions of the Medens instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear.
- Before each procedure, make sure that the driver is properly fitted into the screw or component, otherwise it may occur the deformation, shearing, or even breakage of the driver, screw or component may occur in the application of torque;
- The use of excessive force may compromise the procedure;
- Always use the sequence of Medens products. The use of instruments and / or prosthetic components from other manufacturers and systems does not ensure the perfect function and exempts Medens from offering any product warranty.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instrument or its handling is inadequate, which can cause damage to the surgical sequence of the implant to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The professional must always evaluate the conditions of the instruments after each use. If the instrument is detected to be damaged, markings erased, deformed or excessive wear, must be disposed of in accordance with current regulations, complying with the environmental requirements and the criterion of potentially hazardous products, taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

The Medens bi-digital hexagon driver must be disposed when it loses its functionality or when the user detects excessive wear on the hexagon geometry.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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