IFU

(Instruction for use)

Anilhas / Sleeves / Anilla

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological beings and therefore, to obtain results, it must have been improved to the original standards according to Medens instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral Maxillofacial Surgeons.

DESCRIPTION

As a sleeves, the Guide to Anastomoses is used in conjunction with the Guide to Nursing Anatomy. The guides have an internal hole where drill guides and installation keys are inserted. As Sleeves for a set of 3 pieces: Regular washer Ø5.0mm for installation of Ø3,5 and Ø4,0mm implants, Ø3,7mm narrow washer for the exclusive installation of Ø3,5mm implants and Fixing washer for installation of the pin fixation of the surgical guide in the oral cavity. The Guide Sleeves are made of titanium alloy in accordance with ASTM F136. They feature cylindrical geometry with through-hole and external restraints for Resin Guide Surgical. Invaders as table below:

Diamond surgical drill bits receive the application of a thin biocompatible carbon coating that imparts dark coloration to the cutting surface. This film has as main purposes:

Reference	Model	Color	Stop Diameter (mm)	Inner Diameter (mm)	Body Diameter (mm)	Altura (mm)
AEG-3740	Estreita	Lilás	5,0	3,7	4,0	4
ARG-5060	Regular	Amarela	7,0	5,0	6,0	4
AFG-1380	Fixador	Natural	3,2	1,3	2,6	8

Guide Sleeves are anodized according to the cores in the table above, except for a surgical guide pin washer.

APPLICATION

The Medens Guide surgical technique is indicated for surgical procedures in the jaw or mandible with a guided instrument sequence, allowing the surgery without flap for the installation of implants in edentulous areas.

Guide Sleeves are attached to the Prototyping Surgical Guide. Sistem with the use of Medens Guide instruments and as support to guide the implant during the guided surgery technique. They should be used as per the following table:

Model	Indication of use
Regular	Installation of implants of diameter: 3,3/3,5 e 4,0/4,1 mm
Narrow	Installation of implants of diameter: 3,5mm in regions with little mesiodistal space.
Fixer	Insertion of the surgical guide securing pin

As a Guide to Sleeves, the policies of fabric production and additive manufacturing (prototyping) are polymeric.

CONTRAINDICATIONS

It is contraindicated to reuse according to Guide Sleeves for the accumulation of nerve damage / wear negatively during the surgical procedure.

HANDLING

During the virtual how Sleeves should be selected according to the implant


planned. Regular and distant Sleeves respect the 9 mm distance from the top of the washer to the implant platform. After the manufacture of the Surgical Guide, they must be fixed according to aerial techniques,

used adhesive glue adhesive. After drilling, perform a visual and up-to-date inspection of no load on the top and inside of the washer. Use the products in the Medens Guide line according to the instructions for Use.

* NOTE: The use of ethyl-cyanoacrylate adhesive is recommended.

SANITATION AND STERILIZATION

The Washer Guide is for single use, not extended and should be sanitized and sterilized, following the following recommendations:

PRESENTATION

As a Guide to Sleeves, they are packaged in medical paper packaging and transparent film type. The instrumentation is a process of decontamination that eliminates any type of dirt, but the process of hygiene, described above, must be followed as soon as it is first used.

PRECAUTIONS

- Inappropriate surgical planning may compromise implant and / or instrument performance resulting in system failure, such as loss or fracture.
- Make sure the correct positioning of the Washer Guide and its attachment to the adhesive glue.
- This product is for single use only.
- Reusing this product may cause:
 - *Changes in the physical, mechanical and chemical characteristics, macro and micro structural, original of the product that may impair its intended functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of the related products.*
- During the surgical procedure make sure not to exert lateral movements with the surgical instruments, avoiding the loss of the fixation of the Washer Guide.
- Do not use if the packaging has been tampered with;
- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the Guide Sleeves in their original packaging;
- Always use the product sequence In fi ls. The use of instrumentation and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any warranty of the product.

ADVERSE EFFECTS

Adverse effects will only occur if the choice of instruments is inadequate, which may result in damage to the implant or component to be installed.

STORAGE CONDITIONS

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

DISPOSAL OF MATERIAL

The operator must always evaluate the condition of the instrumentation after each use and it must be discarded in accordance with current regulations, complying with the environmental requirements and the criterion of potentially dangerous products, and if the instrument is damaged, markings erased, deformed and worn. , taking into account the different levels of contamination, physical risks and selective collection.

LIFE CYCLE

Sleeves are for single use and should be discarded shortly after use.

SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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