IFU

(Instruction for use)

Análogo Digital / Digital Analogue / Análogo Digital

PORTUGUÊS

ENGLISH

ESPAÑOL

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

Favor contactar seu distribuidor autorizado.


This product has been designed to incorporate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components in accordance with Emfils instructions and recommendations. All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral and Maxillofacial Surgeons.

DESCRIPTION

The Emfils Digital Analog is part of the family of laboratory components and is used as a replica of the implant that is installed on the patient to perform the laboratory work. In its cervical portion has interface similar to implants, prosthetic components and intermediates. In its apical portion, it has geometry for retention in the prototype (printed) model composed of a Grip system to check "click" fit, a strap for precision seating in the longitudinal axis and parallel to prevent rotation in the cavity. It is made of titanium alloy, according to ASTM F136 standard.

Stop
Greater security in longitudinal fit.

Click
The grip system promotes the interference necessary for the click-fit system.

Parallels
Anti-rotational system in the positioning of the prototype model cavity.

The analogues are presented according to the table below:

Sales Code	Tipo Conexão					Diameter (mm)	Thread	Handle
	CM/LZ	HE	HI	Connection				
				Multi-Unit	Mini 4,8			
ADM-3510	X					Ø3,5	M1,6	MAD-1616
ADM-4010	X					Ø4,0	M1,6	MAD-1616
ADM-5010	X					Ø5,0	M1,6	MAD-1616
ADE-3510		X				Ø3,5	M2,0	MAD-2020
ADE-4110		X				Ø4,1	M2,0	MAD-2020
ADE-5010		X				Ø5,0	M2,0	MAD-2020
ADI-3310			X			Ø3,3	M2,0	MAD-2020
ADI-4010			X			Ø4,0	M2,0	MAD-2020
ADI-5010			X			Ø5,0	M2,0	MAD-2020
ADNR-4010				X		Ø4,0	M2,0	MAD-2020
ADNR-5010				X		Ø5,0	M2,0	MAD-2020
ADMP-4810					X	Ø4,8	M1,4	MAD-1414

APPLICATION

Digital analogs are suitable for simulating dental implant connection and position in a working model in the laboratory phase. Prior to any laboratory manipulation, make sure that the elements are properly secured to the prototyped (printed) model.

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Do not use the digital analogue in plaster models.


Digital analogs are intended for use in prototyped (printed) models only.

Do not reuse the digital analog.

HANDLING

1. Verify CAD/CAM software compatibility and availability at <https://www.emfils.com.br/area-digital/>
2. Make sure the library is up to date as per the instruction attached to the file available for download from the site
3. Obtaining the scan or digital file;

4. Project development in CAD/CAM software compatible with Emfils system;
5. 3d printing of the model;
6. Fit the corresponding knob on the digital analogue (table 1);


7. Bring the digital analog into the prototype model cavity and rotate until it clicks into place (single position).

8. Gently press the digital analogue until it clicks.

9. Inspect the bottom of the prototype model for the position of the digital analog. In addition to the "click" there are other safety devices such as the stop belt and the underside that must be level with the print boundaries.

PRESENTATION

The digital analog is individually wrapped in tyvek® surgical paper packaging and transparent film display. They are supplied ready to use as they undergo a decontamination process that removes any dirt.

STERILIZATION

Supplied non-sterile. Suitable for laboratory use.

PRECAUTIONS

- PRECAUTIONS
- This product is for single use only;
- Reuse may compromise the efficiency of the prosthetic fitting and grip system;
- Digital Analog should be used only with the Emfils digital library compatible with CadCam Software;
- Keep the library updated;
- The Digital Analog should never be altered or modified;
- Do not use if packaging is broken;
- Check the correspondence between the diameter of the analogue and the diameter of the implant installed;
- Careful handling, maintaining the integrity of the platform, avoiding damage and foreign bodies in the plug and screw interface;
- After making the prosthesis, evaluate the appropriate adaptation on the analogue;
- Higher than recommended insertion torques or use of improper connection wrenches can cause material damage and render the system unusable;
- Always use the Emfils product sequence. The use of instruments and / or components from other manufacturers does not guarantee the perfect function of the Emfils implant system and disclaims any product warranty;
- The component must be perfectly seated to the analog before applying the required torque to the screw;
- The digital analogue must be handled by qualified professionals;
- The component must be perfectly seated to the analog before applying the required torque to the screw;
- The digital analogue must be handled by qualified professionals;

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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