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
One of the factors for success in installing a Novo Colosso implant is the use of precise and well-maintained instruments. Medens supports the professional with all the quality and technology that the instrument is manufactured using high quality raw materials and with the utmost precision. Well-maintained instruments do not just mean safety to the patient and the entire surgical team against infection, but it is also a prerequisite for successful treatment! In this topic, we will discuss ways to conserve your instruments, as the maintenance is simple and it is your responsibility to guarantee the quality and durability of all the products developed by Medens.

The instrument cleaning process can be manual or automatic, and involves at least six basic steps: **PRE-CLEANING, DECONTAMINATION, WASHING, RINSING, LUBRICATION AND DRYING.**

1. STAGE: PREPARED CLEANING For manual cleaning, the instrument should be opened, disassembled and dipped, where appropriate, into a suitable container containing water and detergent, preferably enzymatic, at room temperature. It should then be thoroughly rinsed in running water, preferably lukewarm.

2. STAGE: DECONTAMINATION Immediately immerse the instrumentation in an appropriate container containing a disinfectant solution in water, at room temperature (chemical disinfection), or in a heated bath (thermochemical disinfection). The immersion time of the instrument depends on both the operating temperature, the dilution and the type of disinfectant used.

3. STAGE: WASHING The parts should be brushed with a soft bristle brush, paying particular attention to joints, serrations and racks. The joints should remain open, being brushed on both sides, the serrations brushed in both directions towards the serra itself and the racks brushed in both directions towards the teeth. After washing, the instruments should be thoroughly rinsed under running water, and the articulated instruments should be opened and closed a few times during rinsing. It is recommended to use heated water for instrument rinsing.

4. STAGE: LUBRICATION (for collapsible parts). After rinsing, ensure that the instrument is free of any dirt or other debris. All instruments must be lubricated using a water-soluble and non-corrosive lubricant suitable for medical applications, with special attention to joints of the articulated instruments.

5. STAGE: LUBRICATION (for detachable parts). After rinsing, it must be ensured that the instrument is free of any dirt or other residues. Every instrument must be lubricated, using water-soluble and non-corrosive lubricant, suitable for medical application, with special attention to the joints of the articulated instrument.

6. STAGE: DRYING It is recommended:

- When using fabric, which is absorbent and soft;
- each component of a collapsible instrument is dry insulated
- the instrument having a lumen has a completely dry interior.

Note: Special attention should be paid to areas of difficult access where organic tissue retention and deposition of secretions or disinfectant solutions may occur.

IMPORTANT: Water quality is a fundamental factor both for the cleaning process and for the conservation of the instruments. As for example: If in the water contain the excessive presence of chlorides, this will serve to induce the corrosion process of the stainless steel.

RECOMMENDATIONS: It is recommended that all instruments be cleaned immediately after the surgical procedure in which they are used. In case of automatic processes, the manufacturer's instructions must be followed strictly.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms (viruses, bacteria, microbes and fungi), whether in vegetative or sporulated form. Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.

- Under no circumstances should steel straws or other abrasive products be used.

- All surgical instruments and prosthetic components must be sterilized prior to use.

- Sterilization of surgical instruments is not replaced by cleaning. Once sterilized, surgical instruments should only be opened in Surgical Centers under aseptic conditions.

- Parts that have rubbers, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

CARE OF STORAGE AND TRANSPORTATION

It is recommended to keep the surgical instruments away from light, heat and moisture. REF.: Technical Standard NBR 14332

CAUTION: DO NOT STICK LABELS AND DO NOT WRITE IN THE KIT

AUTOCLAVABLE KIT CLEANING INSTRUCTIONS
(Manufacturer's recommendation)

INDICATIONS:

Accommodate instruments for sterilization, handling and transport of the same.

INSTRUCTION FOR USE:

To open the Sterilization Case, the side latches should be disabled with the help of the thumb and forefinger of each hand simultaneously; Open and accommodate the materials you wish to autoclave.

CLEANING:

The removal of organic matter should be done in ultrasonic vats using enzymatic detergent (minimum three enzymes) at dilution recommended by the manufacturer. After washing in running water and observe the presence of residues, repeat the process, if it persists, remove by means of manual cleaning.

DISINFECTION:

Use 0.2% peracetic acid solution for ten minutes.

STERILIZATION:

It should be done by autoclaving, with a temperature of 121 ° C to 135 ° C, observing the recommendations of the autoclave manufacturer's instruction manual. Repeat the process with each use. Caution: Avoid contact with chemicals (solvent, chlorinated hydrocarbons) as this may result in deterioration.

PACKAGING:

Dispose of the product in an exclusive place, in closed cupboards, protected from dust and insects.

COMPOSITION:

Polymer and non-toxic silicone.

PRECAUTION:

Sterilize before use. Do not touch the product on the inner walls of the autoclave.

EXPIRATION DATE:

5 years from the date of manufacture, in unopened packaging. The shelf life of this product depends on the care taken during its use. Successive uses and sterilizations cause material fatigue and limit the life of the material. Check the condition of this product before each use and purchase a new one if you observe excessive fatigue or risk of breakage.